


APPLIX] SMART

It all adds up



ENTERAL NUTRITION

+ Reassurance

- ● ● Easy to operate
- ● ● Splash proof design facilitates an easy cleaning of the pump
- ● ● Connection possibility to nurse call system enabling you be informed immediately
- ● ● Optional keypad lock out prevents unwanted adjustment changes
- ● ● Free-Flow protection with patented tube clamp for added reassurance

+ Flexibility

- ● ● High feeding rates 1 - 600 ml/h
- ● ● Adjustable alarm volume
- ● ● Automatic priming function
- ● ● No pump segment necessary
- ● ● NiMH-Accumulator (reduced memory effect)
- ● ● High mobility (battery life: 24h with 125 ml/h)

+ Comprehensive alarm functions

● ● ● **IT ALL ADDS UP**

Enteral feeding pump for stationary
and ambulatory use

FEEDING PUMP

APPLIX] SMART

the clever, intelligent feeding pump



Clever:
180° rotation of the holder clamp possible



Intelligent:
Patented tube clamp is closed automatically when opening the pump door

ORDERING INFORMATION

	Art.No.
APPLIX SMART	country specific
APPLIX Pump Set Bag	7751711
APPLIX Pump Set Bottle	7751721
APPLIX Pump Set EasyBag	7751731
APPLIX Pump Set VarioLine	7751691
APPLIX Pump Set DuoLine	7752041
APPLIX Infusion Stand	2877851
APPLIX Table Top Pump Stand universal	7751081
APPLIX Adapter for Table Top Pump Stand 7751121	7982101
APPLIX Pole Stand	7982071
APPLIX EasyBag Rucksack	7751361
APPLIX Waistbelt	7750501
APPLIX DuoBag	7902011

7347501 (07/02)