

MEDICAL DEVICES



Base A Data Concentrator

The Base A is a Data Concentrator and supplies mains power for up to 8 Module DPS and/or Module DPS Visio and/or Module MVP.

The Base A enables the connection of Orchestra® Infusion Workstation to a Patient Data Management System via a unique RS 232 cable.



Caring for Life

BASE A

According to the number of intravenous infusions needed, at any combination, up to 8 Modules (Module DPS or Module DPS Visio or Module MVP) can be fitted onto the Base A Orchestra®.

A single mains cord and a single communication cable help reduce installation constraints.

BASE A - TECHNICAL SPECIFICATIONS

MAINS SUPPLY	From 95 to 240 V - 50/60Hz.
BATTERY	Characteristics: 6V 1.1/1.3 Ah - sealed lead acid battery rechargeable. Battery life: 2 hours minimum.
ALARMS, PREALARMS AND SAFETY	The Base A gives a synthetic view of the operation of connected Modules.
COMPLIANCE, STANDARDS	Complies with the 93/42 EEC European Medical Device Directive: CE 0459 Complies with EN 60601-1, EN 60601-1-2 and EN 60601-2-24 requirements.
WATERPROOFNESS	Protection against ingress of liquid IP34 (CEI 529).
ELECTRICAL SAFETY	Protection against leakage current: type CF equipment. Protection against electric shocks: class II equipment.
COMMUNICATION INTERFACE	Compatible with major PDMS software.
RS232 INTERFACE	One unique cable for all connected Modules.
DIMENSIONS / WEIGHT	W=55, H=315, D=200 mm / 2 kg.

Orchestra®

Infusion Workstation composed of Bases and Infusion Modules adapted to every clinical practice.

Module DPS
Modular Syringe Pump.



Module DPS Visio
Modular Syringe Pump.



Module MVP
Modular Volumetric Pump.



Base A
Data concentrator.



Base Intensive
Unique 2 channel relay function, Customised drug library, Infusion history management.



Base Primea®
Enables simultaneous control of both hypnotic and analgesic effects in TCI.



Orchestra® Mobile stand

Whatever may be the architecture of your department, Orchestra® Mobile Stand provides the choice of customised accessories.



Disposable

A wide range for use with Module DPS and Module MVP. Catalogue is available from our Customer Service.

Please Note: Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please consult us for the current situation.